


THE BIPOLAR 5 mm OFFICE RESECTOSCOPE



Click or scan here
to watch the teaser video

Versatile use









- Minimal outer diameter
 - Possible to perform resection in the doctor's office or clinical outpatient setting, even without anesthesia
 - Little or even no dilation required (depending on the patient)
- Optional use as „see & treat“ hysteroscope
- Real bipolar system with current not returned via the sheath




STORZ
KARL STORZ – ENDOSKOPE


THE BIPOLAR 5 mm OFFICE RESECTOSCOPE

Due to its small diameter, this resectoscope offers the possibility to perform resections outside the OR. Infertility patients or those suffering from abnormal bleeding can thus be treated in the doctor's office or clinical outpatient setting, even without anesthesia.

Optic and Working Element Set	
26120AA	HOPKINS® Straight Forward Telescope 0° , diameter 2.9 mm, length 30 cm
26053SCK	Resectoscope Sheath , 15 Fr., oblique beak, rotating inner sheath, color code: green
26053OC	Standard Obturator
26053EBH	Working Element Set , bipolar
26053CD	Telescope Bridge , with channel for semi rigid 5 Fr. operating instruments, for use with Resectoscope Sheath 26053SC

HF Electrodes for various indications (sterile, for single use, package of 10)			
	011050-10	Cutting Loop , bipolar, 15 Fr.,	STERILE 
	011051-10	Cutting Electrode , bipolar, pointed, 15 Fr.,	STERILE 
	011052-10	Coagulation Electrode , bipolar, ball end, 15 Fr.,	STERILE 
	011053-10	Cutting Loop , bipolar, longitudinal, 15 Fr.,	STERILE 

Cold Loops for blunt dissection/enucleation or retrieving resected tissue		
	26053R	MAZZON Cold Loop , straight, 15 Fr.
	26053T	MAZZON Cold Loop , rake-shaped, 15 Fr.
	26053U	MAZZON Cold Loop , knife-shaped, 15 Fr.

HF Device		Pump System	
UH400	AUTOCON® III 400 , high-end HF generator, without CE certification, power supply 220--240 VAC, 50/60 Hz	UP210	ENDOMAT® SELECT
		UP602	HYSTEROSCOPY Software
UH400E	Same , with CE certification	031523-10	Tubing Set for Irrigation , pressure-controlled STERILE 

For direct contact email us at gyn@karlstorz.com

Es wird empfohlen, vor der Verwendung die Eignung der Produkte für den geplanten Eingriff zu überprüfen. Bitte beachten, dass die in diesem Medium aufgeführten Produkte möglicherweise aufgrund unterschiedlicher Zulassungsvoraussetzungen noch nicht in jedem Land verfügbar sind.

STORZ
KARL STORZ — ENDOSKOPE
THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com